


# CaseCenter Login

User name

Password



**3DHISTECH**

*Solution resides in the details*

## CASECENTER


**GET FULL CONTROL ON YOUR  
DIGITAL SLIDES!**

CaseCenter is a full featured digital slide management software from 3DHISTECH. Its flexible structure can be adapted to various fields, including clinical pathology, research applications, teleconsultation and education. Integration with existing medical information systems is also possible.



## DIGITAL SLIDE MANAGEMENT

- Flexible folder and case structure
- Use barcodes to organize your digital slides, macro images and project files easily
- Multiple user levels for different access to information
- Storage can grow with your needs without limitation on the number of slides stored
- Free API for custom solutions



## APPLICATIONS

### Clinical & Research


- Education
- Consultation
- Second opinion

Integration with HIS/LIS systems available

Faster diagnosis by automatic alignment of different staining of the same block


### VirtualTray

- All slides of a case at once in full resolution
- Slides or slide parts can be replicated without changing the original slide or increasing storage needs



### Teleconsultation

- Over 100 simultaneous viewers possible
- Session setup in less than a minute
- Session notification with in every common calendar software
- Share control of the slide
- Built-in chat function
- Slides are available before and after the Teleconsultation session for review



### Slide archiving

Secure database with automated backup and integrity check

Production and development by

**3DHISTECH LTD.**  
Konkoly-Thege Miklos st. 29-33., Bldg. 6.  
H-1121 Budapest, Hungary  
Phone: +36-1-392-2274 Fax: +36-1-392-2723  
Info: info@3dhistech.com www.3dhistech.com



**3DHISTECH**