



## QuantCenter 2.0 Main features




### IHC quantification

Applying **PatternQuant** presegmentation you can detect the cancer region accurately. Background or connective tissue region can be eliminated and you can focus only on the region of interest and run cell based measurements only on the cancer region.




### Measurement Visualization

By using the Measurement Visualization function you can filter and display measurement results either in Piechart, Scatterplot, Histogram, or in the Gallery.



### Review functions

The Gallery view allows you to apply a correct relocalization on digital slides. You can modify object detection by using the edit shape or the rescore function.



### Export

Beside the sophisticated export to Excel, all of the Gallery or diagram images can be exported to a variety of common final formats. Diagnostic modules are supported with a printable measurement report function.