


CaseCenter 2.4




- **Faster interface**
Manage your cases and digital slides smoothly

- **Automatic slide alignment**
Automatic and rapid alignment of series of sections even with different stainings. It is not necessary to rotate digital slides manually.



- **VirtualTray**
All slides of a case represented in one view by one single click for faster case reading and slide comparison. Combine multiple digital slides or regions of interest into one VirtualTray without the need for more storage!



- **TMA workflow support**
from target marking up to the core visualization. Use CaseCenter's remote slide access for sample designation or TMA slide evaluation: fast, convenient, and easy.

- **Transparent login method for high traffic user authentication**

Login queue and automatic logout. It is possible for the administrator to log users out after a specified time of inactivity. A login queue handles the shortage of free licenses. Users get in the queue according to their user levels (users with higher rank are put forward in the queue) and their queue number (position in the queue) is displayed and is automatically updated.